

HC100™ Z-Band® Wristbands



Zebra's Z-Band® direct thermal wristbands allow you to print patient information and bar codes directly on the wristbands, so there is no need for labels and the time-consuming wristband assembly inherent to laser technology.

In addition, Zebra's Z-Band direct thermal wristbands:

- Feature a silver antimicrobial coating that protects the wristband.
- Are scannable for up to 14 days, promoting quality patient care and safety.
- Are quick and easy for staff to print and securely fasten in less than a minute using an adhesive tab or color clip.
- Were tested by an independent MRI testing company and were determined **"MR-Safe" (MR Safe—an item that poses no known**








hazards in all MRI environments) according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08: Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. These wristbands do not appear to present an additional hazard or risk to a patient undergoing an MRI procedure.



Zebra Z-Band wristbands are available with adhesive or clip closure, in a variety of sizes and colors to meet your organization's needs.


Z-Band Comfort

- Adhesive closure.
- Soft and flexible.
- Superior resistance to hand sanitizers such as Purell®.
- 50 percent more wristbands per cartridge than Z-Band Direct.

Part Number	Size	Color	
Adult			
10010951K	1" x 11"	White	
10010951-1K (Red), 10010951-2K (Blue), 10010951-3K (Green), 10010951-4K (Purple), 10010951-5K (Pink), 10010951-6K (Orange), 10010951-7K (Yellow)			
10011952K	0.75" x 11"	White	
10011963K	1.1875" x 11"	White	
Pediatric			
10011954K	1" x 7"	White	
Infant			
10011953K	1" x 6"	White	



Soft Infant Z-Band Direct

- Adhesive closure.
- Soft nylon backing is perfect for delicate skin.
- Adjustable band fits the smallest patients.
- Flat print area ensures ease of scanning.

10007746K	1" x 7.6875"	White	
-----------	--------------	-------	---

Z-Band QuickClip™

- Clip closure.
- Color clips available in red, yellow, purple, green, pink and blue.

Adult			
10007000K	1.1875" x 11"	White	
Infant			
10007003K	1" x 7"	White	

A bar code that fails to scan or requires repeated attempts to scan jeopardizes patient safety and staff productivity. Zebra's supplies research and development team tested the image durability of Z-Band direct thermal wristbands to ensure optimum image durability and scanning performance when subjected to common solvents found in hospitals including 91% Isopropyl Alcohol, 70% Ethyl Rubbing Alcohol, Betadine and hand sanitizer, as well as abrasion and water.

For more details on the durability laboratory testing performed, please download our white paper on www.zebra.com titled "Laboratory Testing of Zebra's Z-Band® Direct Thermal Antimicrobial Wristbands."

©2010 ZIH Corp. All product names and numbers are Zebra trademarks, and Zebra, the Zebra head graphic, and Z-Band are all registered trademarks of ZIH Corp. All rights reserved. Purell is a registered trademark of Warner-Lambert Company. All other trademarks are property of their respective owners.



www.zebra.com

Corporate Headquarters
+1 800 423 0442
E-mail: inquiry4@zebra.com

Asia-Pacific Headquarters
+65 6858 0722
E-mail: apacchannelmarketing@zebra.com

EMEA Headquarters
+44 (0)1628 556000
E-mail: ms europe@zebra.com

Latin America Headquarters
+1 847 955 2283
E-mail: inquiry4@zebra.com

Locations

USA: California, Georgia, Rhode Island, Texas, Wisconsin **Europe:** France, Germany, Italy, Netherlands, Poland, Spain, Sweden **Asia Pacific:** Australia, China, India, Japan, South Korea **Latin America:** Argentina, Brazil, Florida (USA), Mexico **Africa/Middle East:** Russia, South Africa, United Arab Emirates

GSA#: GS-35F-0268N
P1015680 Rev. 3 (8/10)





www.L-TronDirect.com

800-830-9523

info@L-Tron.com

596 Fishers Station Dr | Victor, NY | 14564 | Suite 1 A

www.L-Tron.com

